



VANDERBILT  
UNIVERSITY

# Porous Silicon Flakes on a Flexible Substrate for Real-Time Biosensing using Smartphone Technology

Caitlin Carfano<sup>1</sup>, Tengfei Cao<sup>2</sup>, and Sharon M. Weiss<sup>3</sup>

<sup>1</sup> Department of Electrical and Computer Engineering, The George Washington University, Washington, DC

<sup>2</sup> Interdisciplinary Graduate Program in Materials Science, Vanderbilt University, Nashville, TN

<sup>3</sup> Department of Electrical Engineering and Computer Science, Vanderbilt University, Nashville, TN




## Motivation

### Low Cost Diagnostics:

- Cannot rely on expensive lab equipment
- Must not take long time to get results
- Especially needed in developing world
- Smartphones are ubiquitous and can be easily implemented in a point-of-care system

### Meeting ASSURED Criteria:

- ✓ Affordable
- ✓ Sensitive
- ✓ Specific
- ✓ User-Friendly
- ✓ Rapid and Robust
- ✓ Equipment Free
- ✓ Deliverable to Users




[1] Figure: More than half the deaths in the poorest countries are a result of infectious and parasitic diseases

## Porous Silicon as a Biosensor

### [2] Porous Silicon (PSi) Properties:



- Large internal surface area ( $> 100 \text{ m}^2/\text{cm}^3$ )
- Widely tunable pore size ( $\sim 1 - 100 \text{ nm}$ )
- Easy to fabricate (electrochemical etching)



### Operating Principle:



- Multilayer PSi film selectively reflects certain colors of light
- When pores are partially or completely filled, the reflected color changes

### Cross-Sectional View





### Porous Silicon Particles:

Single  $\sim 30 \mu\text{m}$  particles broken from a PSi film show a reflectance spectrum that directly correlates to the color perceived by the eye





## Device Fabrication


1 Attach porous silicon to a flat layer of PDMS; place tape on filter paper and spread PDMS across it



2 Peel tape off the paper and place the porous silicon/PDMS unit to paper/PDMS. Bake at high temperature in oven




3 Place device paper first in solution and observe color change




## Future Work

- Change fabrication steps to prevent hydrophobic PDMS from entering the pores of the porous silicon





- Create a user-friendly app to test different analytes simultaneously

Proposed app design



## Sensing with Glucose



## Conclusions

- Demonstrated that PSi can be transferred directly onto flexible PDMS or plastic films without compromising spectral integrity
- Presented new approach for low-cost, easy-to-use biosensor using a combination of PDMS, paper, and hydrophobic barriers with PSi and smartphone
- Demonstrated that using RGB data from a smartphone camera can be used instead of a traditional spectrometer to accurately measure sensing events in PSi

## References:

- [1] Yager, Paul, et al. "Microfluidic diagnostic technologies for global public health." *Nature* 442.7101 (2006): 412.
- [2] Rong, G., et al. "Label-free porous silicon membrane waveguide for DNA sensing." *Applied Physics Letters* 93.16 (2008): 161109.
- [3] Parolo, C. and Arben, M. "based nanobiosensors for diagnostics." *Chemical Society Reviews* 42.2 (2013): 450-457.
- [4] Costa, M. N., et al. "A low cost, safe, disposable, rapid and self-sustainable paper-based platform for diagnostic testing: lab-on-paper." *Nanotechnology* 25.9 (2014): 094006.



**Acknowledgements:** This research was supported in part by the NSF ViNSE REU grant DMR-1560414 and the ARO grant W911NF-15-1-0176. I gratefully acknowledge Alice Leach for helping in the clean room, Tianjiao Wang for single particle reflectance measurements in Prof. Yaqiong Xu's laboratory, as well as everybody in the Weiss lab for support.

