



Porous Silicon on Glass for Low-Cost Diagnostics

Crystal A. Nattoo¹, Tengfei Cao², and Sharon M. Weiss^{2,3}



1. Electrical and Computer Engineering Department, University of Miami, Coral Gables, FL
2. Interdisciplinary Graduate Program in Materials Science, Vanderbilt University, Nashville, TN
3. Department of Electrical Engineering and Computer Science, Vanderbilt University, Nashville, TN

Background

Current System of Diagnostic Testing: Chips-In-A-Lab

- Use large and expensive equipment
- Samples usually processed far from the patient
- Difficult to track down patient for treatment in rural areas


Characteristics of the ideal diagnostic test — ASSURED

- **Affordable** by those at risk of infection.
- **Sensitive** (few false-negatives).
- **Specific** (few false-positives).
- **User-friendly** (simple to perform and requiring minimal training).
- **Rapid** (to enable treatment at first visit) and **Robust** (does not require refrigerated storage).
- **Equipment-free**.
- **Delivered** to those who need it.

Reference: [1]

New System: Lab-On-A-Chip


- Cost-effective, fast, and portable
- Use a smartphone for the computing, user interface, and sensing



Materials


Porous Silicon (PSi) Film

- **Microcavity filter:** thickness 1.5 - 3 μm
- **Pore diameter:** 20 - 30 nm
- **Fabrication:** electrochemical etch and liftoff followed by chemical oxidation
- Large surface area $\approx 100 \text{ m}^2/\text{cm}^3$
- Easy surface modification





Glass

- Enables robust platform



Paper

- **Filter paper**
 - Facilitates molecule transportations
 - Blocks undesired molecules
- **Blotting paper**
 - Works as a pump to draw liquid across filter paper





Progression of Film Quality through Fabrication Process

Testing With Biomolecules


- Measure reflectance spectra using Ocean Optics® spectrometer
- Reflectance spectrum red-shifts when analyte infiltrates into the pores
- Magnitude of shift quantifies the molecules added




Flow Cell




Biomolecule Tests w/ Concentrations of Glucose





Glucose Test (Zoomed)



Proof of Reversibility





Reversibility (Zoomed)




Sensing with Smartphone Camera

- Use smartphone to take videos of the sample
- Film color changes when analyte infiltrates into the pores
- Magnitude of change in light intensity at the red channel (R value) quantifies the molecules added
- Film spectrum shifts from around 500 nm (green) to 600 nm (red-orange) when wet.



Contour map of R values obtained from video of PSi film during experiment




Conclusion

- Successfully lifted porous silicon thin film off silicon wafer and transferred onto glass slide without compromising optical properties
- Built a robust PSi sensing platform
- Measured change in porous silicon thin film refractive index due to infiltration with glucose using both a spectrometer and a smartphone camera


Future Work

- Correlate spectral shift to smartphone readings
- Test new platform designs to reduce inherent error due to motion

Current Platform



New Platform



References

- [1] Mabey, David, et al. "Diagnostics for the developing world." *Nature reviews. Microbiology* 2.3 (2004): 231.
- [2] Erickson, David, et al. "Smartphone technology can be transformative to the deployment of lab-on-chip diagnostics." *Lab on a Chip* 14.17 (2014): 3159-3164.
- [3] Sailor, Michael J. *Porous silicon in practice: preparation, characterization and applications*. John Wiley & Sons, 2012.

Acknowledgments

I would like to acknowledge the Weiss Lab for providing helpful feedback while I worked on this project.

This work was supported in part by the NSF REU grant DMR-1560414 and the Army Research Office grant W911NF-15-1-0176.