

Kaleigh A. Ellis ¹, Kemar R. Reid ², Sandra J. Rosenthal ²

¹ Department of Chemistry and Physics, Saint Mary's College, Notre Dame, IN

² Vanderbilt Institute of Nanoscale Science and Engineering, Vanderbilt University, Nashville, TN

Background

What are Quantum Dots?

- Semiconductor nanocrystals (2-15 nm) with size-tunable optical properties
- Larger QDs have a narrower band gap, resulting in lower-energy photon emissions ("red shift")
- Applications include solid state lighting, solar cells, and **single-protein imaging and tracking**




Figure 1: Quantum confinement in nanocrystals of different sizes ¹

Core/Shell Quantum Dots

- Shelling QD cores with inorganic material (e.g. CdS) confines the **electron-hole pair** inside the QD
- Shelled cores have increased photostability and quantum yield




Figure 2: Quantum confinement of electron (e⁻) and hole (h⁺) wavefunctions between conduction (CB) and valence (VB) bands in a core/shell system

Ideal QDs for Biological Imaging

- Photostable**
 - "Non-blinking"
- Bright**
- High quantum yield
- Narrow Emission Peak**
- Small size-distribution

Figure 3: Laser-scanning confocal microscopy of a HEK293 cell (scale bar 10 μm) ²
Yellow: Human dopamine transporter fused to yellow fluorescent protein at the N-terminus
Red: QDs (655 nm max emission)

Research Goals

Engineer stable, bright, and size-uniform QDs for single-protein imaging by:

- Growing a CdS shell onto alloyed CdSSe cores




Figure 4: Reaction scheme for synthesis of CdSe-CdS core alloy ³

- Studying photoluminescence and structural properties
- Optimizing the quantum yield and size-distribution:**
 - Quantum yield is determined via fluorescence spectroscopy
 - Particle size obtained from transmission electron microscopy (TEM) images
- Two shelling procedures were performed and compared

Methods and Results

Successive Ion Layer Adsorption and Reaction (SILAR)




Figure 5: Reaction conditions for SILAR method




Figure 8: Normalized emission spectra for SILAR method




Figure 6: TEM images for cores and SILAR-shelled QDs (scale bars 50 nm) ⁴




Figure 7: Red-shift of SILAR-shelled QDs under UV-light. From left to right: Cores, 1 ML, 3 ML, 5 ML

Table 1: Summary of results for SILAR method

Sample	Quantum Yield (%)	Max Emission (nm)	FWHM (nm)
Cores	74.3	579.16	39.58
1 ML	54.3	586.56	41.90
3 ML	17.4	611.38	40.95
5 ML	33.0	621.88	38.08

Continuous Precursor Infusion




Figure 9: Reaction conditions for continuous infusion method




Figure 12: Normalized emission spectra for continuous infusion method




Figure 10: TEM images for cores and continuous infusion-shelled QDs (scale bars 10 nm) ⁴




Figure 11: Red-shift of continuous infusion-shelled QDs under UV-light. From left to right: Cores, 40 min, 75 min, 90 min

Table 2: Summary of results for continuous infusion method

Sample	Quantum Yield (%)	Max Emission (nm)	FWHM (nm)
Cores	74.3	589.59	40.09
40 min	54.0	608.09	34.97
75 min	45.9	615.50	34.53
90 min	62.0	615.15	34.51

Conclusions

SILAR Method

- Uneven size distribution
- Lower quantum yield
- Long reaction time

Continuous Infusion Method

- Narrower size distribution
- High quantum yield
- Shorter reaction time

Applications in Biological Imaging

- A thicker CdS shell:
 - Makes the QD more robust to biological ligand exchanges and conjugation
 - Protects the QD from a harsh biological environment
 - Confines the electron-hole pair inside the QD and retains optical properties




Figure 13: Core/Shell QD under ligand-exchange and biological conjugation for single-protein tracking ⁵

Future Work

Further optimize shelling procedure

- Use syringe pumps for more consistent precursor injection
- Size-separate alloyed cores via column chromatography to produce better shelled-QD size distribution
- Study other optical properties

Photostability

- "Blinking" and correlation to QD structure




Figure 14: Correlation between photostability and QD structure ⁶
a) Full fluorescence intensity
b) Corresponding histograms
c) TEM of individual QDs (scale bars 2 nm)

References

- Image from <http://www.sigmaaldrich.com/content/dam/sigma-aldrich/materials-science/nanomaterials/quantum-confinement-effect.jpg>
- Image courtesy of Dr. Oleg Kovtun
- M.A. Harrison, *Isr. J. Chem.* **2012**, *52*, 1063-1072
- Images obtained by Kemar Reid
- J.C. Chang, *J. Neurosci.* **2012**, *32*(26), 8919-8929
- N.J. Orfield, *ACS Nano.* **2015**, *9*, 831-839
- SMC logo obtained from https://www.saintmarys.edu/files/images/SMCsq_Logo_rgb_0.jpg
- VU logo obtained from <http://news.vanderbilt.edu/files/vu06b.jpg>
- NSF logo obtained from <http://www.nsf.gov/images/logos/nsf1.jpg>
- VINSE logo obtained from <http://www.vanderbilt.edu/vinse>

Acknowledgements

Special thanks to the National Science Foundation (grant DMR-1263182) and the Vanderbilt Institute of Nanoscale Science and Engineering (VINSE) for this REU opportunity. Additional thanks to all of the members of the Rosenthal research lab (Dr. Sandra Rosenthal, Dr. James McBride, Dr. Ian Tomlinson, Dr. Oleg Kovtun, Nathaniel Freymeyer, Talitha Frecker, Louie Thal, Cayla Ontko, and especially Kemar Reid) for all their assistance and guidance. I am also grateful for the faculty and staff of the Saint Mary's College Department of Chemistry and Physics for all of their instruction and mentorship, especially Drs. Toni Barstis and Chris Dunlap for their assistance in applying for this program.

